


BOTTOM VIEW

ALL DIMENSIONS ARE IN INCHES

1/4-20 THREAD THRU
A .33 THICK PLATE

M6 THREAD THRU
A .33 THICK PLATE,
3 PLACES

1/4-20 THREAD THRU
A .33 THICK PLATE,
3 PLACES



X AXIS LOCK
FITS 5/64 OR 2.0mm HEX

Y AXIS DIFFERENTIAL SCREW
ADJUSTMENT, FITS 5/64 OR 2.0mm HEX

Y AXIS ADJUSTMENT
FITS 5/64 OR
2.0mm HEX

Y AXIS COARSE ADJUSTMENT LOCK
FITS 5/64 OR 2.0mm HEX

Y AXIS LOCK
FITS 5/64 OR 2.0mm HEX

X AXIS COARSE ADJUSTMENT LOCK
FITS 5/64 OR 2.0mm HEX

X AXIS DIFFERENTIAL SCREW
ADJUSTMENT, FITS 5/64 OR 2.0mm HEX

CLEARANCE HOLE FOR
1/4-20 OR M6 CAP SCREW
2 PLACES

X AXIS ADJUSTMENT
FITS 5/64 OR 2.0mm HEX

1.50 OPTIC
AXIS HEIGHT TYP.

Ø.681
THRU

FOCUS ADJUSTMENT LOCK
FITS 5/64 OR 2.0mm HEX,
6 PLACES

